

Datasheet: 4C002

Description:	FOUR-COLOR HUMAN CD45RA/CD45RO/CD3/CD4 FLOW KIT
Specificity:	CD45RA/CD45RO/CD3/CD4
Format:	4 Color
Product Type:	Four Color Reagent
Clone:	F8-11-13 / UCHL1 / UCHT1 / RPA-T4
Isotype:	Cocktail
Quantity:	50 TESTS/0.5ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Product Form

Four color combination consisting of FITC, RPE, RPE-Cy5.5 and APC conjugated monoclonal antibodies mixed in optimal ratio - lyophilised.

FITC conjugated CD45RA (Mouse IgG1)
 RPE conjugated CD45RO (Mouse IgG2a)
 RPE-Cy5.5 conjugated CD3 (Mouse IgG1)
 APC conjugated CD4 (Mouse IgG1)

Reconstitution

Reconstitute with 0.5ml distilled water

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	APC	650	661
	FITC	490	525
	RPE 488nm laser	496	578
	RPE 561nm laser	546	578
	RPE-Cy5.5 488nm laser	496	695
	RPE-Cy5.5 561nm laser	546	695

Buffer Solution

Phosphate buffered saline

**Preservative
Stabilisers**

0.09% Sodium Azide (NaN₃)
 1% Bovine Serum Albumin
 5% Sucrose

External Database Links	UniProt: P08575 Related reagents P07766 Related reagents P01730 Related reagents
	Entrez Gene: 5788 PTPRC Related reagents 916 CD3E Related reagents 920 CD4 Related reagents
Synonyms	CD45, T3E
Specificity	<p>Four-Color Human CD45RA/CD45RO/CD3/CD4 Flow Kit can be used for single-step identification of human helper/inducer naïve (CD3+CD4+CD45RA+) and helper/inducer memory (CD3+CD4+CD45RO+) T-cell subsets.</p> <p>CD3 is a member of the immunoglobulin superfamily, which acts as a mediator of signal transduction, through association with the α/β or γ/δ T-cell receptor (TCR). Mammalian CD3 is a multimeric protein composed of four distinct polypeptide chains (ε, γ, δ, ζ), consisting of two heterodimers (εγ, εδ) and one homodimer (ζζ). CD3 is expressed by a high-percentage of circulating peripheral T-cells and is considered a pan T-cell marker. Clone UCHT1 specifically recognizes the ~20 kDa CD3ε chain.</p> <p>CD4 is a single chain transmembrane glycoprotein and member of the immunoglobulin superfamily, which interacts directly with class II molecules of the major histocompatibility complex (MHCII), and is the primary receptor for the human immunodeficiency virus (HIV). CD4 is expressed by the helper/inducer (CD3+CD4+) T-cell subset, the majority of thymocytes, and at low levels by monocytes and tissue macrophages. The CD4 lymphocyte count provides an excellent indication of immune system health. Clone RPA-T4 binds to the D1 domain (CDR1 and CDR3 epitopes) of the CD4 antigen.</p> <p>CD45RA is an isoform of CD45, a member of the protein tyrosine phosphatase family, otherwise known as leucocyte common antigen (LCA), which acts as a positive regulator of T-cell co-activation, and is essential for the activation of T-cells through the antigen receptor. CD45RA is expressed by naïve T-cells, but not by activated and memory T-cells. Consequently, clone F8-11-13 can be used to distinguish between these cell populations.</p> <p>CD45RO is the low molecular weight isoform (~180 kDa) of CD45, the leucocyte common antigen (LCA). CD45RO is expressed by approximately 40% of resting peripheral blood T-cells, the majority of activated and memory T-cells and most thymocytes. It is also expressed on macrophages, monocytes and granulocytes. In combination with clone F8-11-13, clone UCHL1 can be used to discern naïve and memory T-cells.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
Storage	Prior to reconstitution store at +4°C. After reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation #10075 available at: 10075: https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

FLOW CYTOMETRY ABSOLUTE COUNT STANDARD™ (FCSC580)

ERYTHROLYSE RED BLOOD CELL LYSING BUFFER (10x) (BUF04B)

ERYTHROLYSE RED BLOOD CELL LYSING BUFFER (10x) (BUF04C)

HUMAN SEROBLOCK (BUF070A)

HUMAN SEROBLOCK (BUF070B)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M308886:170816'

Printed on 02 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)